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- Remove all pads and straps from the metal frame. Close all velcro straps into material prior to laundering. Clean metal frame with lukewarm water and mild soap.
- Do NOT use this brace during high activity levels (active sports, running etc.).

- Extra padding is available as an accessory. Contact your fitter for further information.
- Do check the condition of the leg frequently.
- Do check the condition of the brace before use every day and contact your orthotist if any abnormal conditions appear.

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GB

Make sure all straps on the back side of the brace are closed and the pad is located in the correct position as shown above. Then place the brace (with knee flexed in 45°) on the back side of the leg. If the brace includes side pads, make sure they are properly placed between the hinge on the brace and your leg.

Start by securing the first strap below the knee (1).

Next, secure the bottom strap (2).

Secure the strap just above the knee (3).

Last, secure the top strap (4).

Check that the brace fits comfortably on the leg.

To remove, open only one side of the front straps. Leave the back straps closed.

If the brace has metal loops on one side this is where the brace should be opened for doffing and donning.

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